 **16 min**
96 Samples

Swift96

Automatic Nucleic Acid Extraction Instrument (Model: DA8045)

CE

NMPA

FSC

Swift96

Swift96 is an automated nucleic extraction instrument adopting magnetic bead technology. It is used for the extraction, purification of nucleic acid from sample in laboratory environment. Together with the nucleic acid extraction kit, Swift96 is able to extract and purify nucleic acid from 96 samples in 16 minutes.

Features


User-friendly

Convenient setting, visualized interface, touch screen and easy to operate.

Motion control

Precise and intelligent motion control can greatly reduce the loss of magnetic beads.

Free programming


A powerful program edit function can meet requirements of different reagents.


Safe

UV sterilization

Structure


Applicable reagent: DA1003 (96 tests/kit)

Operation Procedure


Magnetic Bead Technology Based \approx 16 minutes / 96 samples (Use DA1003 Purification Kit)

1


Select the program and set up

2

Place reagent plates, magnetic bar sleeve and consumables

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NMPA

中华人民共和国
药品GMP证书
 CERTIFICATE OF GOOD MANUFACTURING PRACTICES FOR PHARMACEUTICAL PRODUCTS
 PEOPLE'S REPUBLIC OF CHINA

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经审查, 符合中华人民共和国《药品生产质量管理规范》要求。
 特此发证。
 This is to certify that the above-mentioned manufacturer complies with the requirements of Chinese Good Manufacturing Practices for Pharmaceutical Products.

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GMP

 By Royal Charter

bsi.
Certificate of Registration
 QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

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Holds Certificate No: **MD 98885**
 and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:
 The design, manufacture and sales of reagents for use in nucleic acid and immunological based in vitro diagnostic testing.
 基于核酸和免疫学体外诊断试剂的设计、制造和销售。

 Gary E Slack, Senior Vice President - Medical Devices

For and on behalf of BSI:
 Original Registration Date: 2006-11-22 Effective Date: 2021-11-04
 Latest Revision Date: 2021-10-29 Expiry Date: 2024-11-03
 Page 1 of 2

...making excellence a habit™

This certificate was issued electronically and retains the integrity of BSI and is bound by the conditions of contract.
 An electronic certificate can be authenticated online.
 Printed copies can be validated at www.bsi.com/CheckCertificate or telephone +44 (0) 800 9000.
 Information and Contact: BSI, 389 Chiswick Court, Uxbridge, Middlesex, UK. Tel: +44 (0) 800 9000
 BSI is a registered charity, registered in England under number 3090712 and in the UK Charity Commission under number 1101204.

ISO13485

Parameters

Throughput	Up to 96 samples
Reaction volume	20-1000 µL
Anti-pollution design	UV light, Exhaust rate > 1.2 m/s
Temperature range	Room temperature to 99°C
Communication interface	Ethernet port
Magnetic beads recovery rate	≥98%
Dimension (L×W×H)	560×580×635 mm
Net weight (kg)	62
Certificates	CE, NMPA, FSC

Daan Gene Co., Ltd.

CHINA (Headquarters)

📍 **Add:** No.19, Xiangshan Road, Science Park, High & New Technology Development District, Guangzhou, Guangdong, P.R. China
 ☎ **Tel.:** +86-20-32068126
 📠 **Fax:** +86-20-32068352
 ✉ **Email:** marketing@daangene.com

CANADA office

📍 **Add.:** #413-4538 Kingsway, Burnaby, BC V5H 4T9, Canada
 ☎ **Tel.:** +16044517588
 📠 **Fax:** +604451.7587
 ✉ **Email:** gbai@daangene.net

COLOMBIA office

📍 **Add.:** Cl 93 No.16 46 Of 702, Bogotá D.C.
 ☎ **Tel.:** +57 3103940417
 ✉ **Email:** alonso.zheng@daangene.net

🌐 **Website:** <https://en.daangene.com/>