




SAMPLE COLLECTION



SAMPLE LOADING AND DETECTION



RESULT



CONTACT US

CHINA (Headquarters)

Address: No.19, Xiangshan Road, Science Park, High & New Technology Development District, Guangzhou, Guangdong, P.R. China
Tel.: +86-20-32068126
Fax: +86-20-32068352
Email: marketing@daangene.com

CANADA office

Address.: #413-4538 Kingsway, Burnaby, BC V5H 4T9, Canada
Tel.: +1604 451 7588
Fax: +604451.7587
Email: gбай@daangene.net

COLOMBIA office

Add.: Cl 93 No.16 46 Of 702, Bogotá D.C.
Tel.: +57 3103940417
Email: alonso.zheng@daangene.net

COVID-19 Antigen Rapid Testing Kit

For home/self-testing use



Daan Gene Co., Ltd.
Website: <https://en.daangene.com/>



FEATURES

- 1.High sensitivity:** 100 TCID₅₀/mL
- 2.High precision:** Detection range covers wild type SARS-CoV-2 and its variants.
- 3.Fast:** Get results within 10-15 mins.
- 4.Easy to operate:** Only 3 steps: sampling, dilution and sample loading.
- 5.Easy to result:** Interpret the fluorescent result by naked eye.

SPECIFICATION

【Packaging】

1 Test/Kit & 10 Tests/Kit

【Specimen type】

Anterior nasal specimens

【Storage condition】

Unopened and stored at 2-30°C

【Shelf life】

12 months

PRODUCT INFORMATION

COVID-19 Antigen Rapid Testing Kit is an immunochromatographic reagent based on the principle of the double-antibody sandwich method. This kit is intended for the in vitro qualitative detection of SARS-CoV-2 antigen in human anterior nasal swab specimens from individuals who are asymptomatic or suspected of SARS-CoV-2 infection.

